

# CE Diameter Report - Volume Distribution



## Measurement Details

**Sample Name:** Bactiblock 101

**SOP Name:** MATERIE PRIME.vsop

**Edited:** True

**User Name:** morfologi

**Edited On:** 22 April 2014 18:41:25

**Date:** 22 April 2014 18:34:02

## Analysis

**SOP Analysis:** 3.0

**SOP Fill Holes:** True

**SOP Segmentation:** None

**SOP Trash Size:** 10

**SOP Optic(s) used:** 10x, 50x

## Result

**CE Diameter Minimum ( $\mu\text{m}$ ):** 0.22

**CE Diameter D[v, 0.1]:** 3.968

**CE Diameter Maximum ( $\mu\text{m}$ ):** 90.60

**CE Diameter D[v, 0.5]:** 12.75

**CE Diameter D[4,3] ( $\mu\text{m}$ ):** 22.88

**CE Diameter D[v, 0.9]:** 89.24


**CE Diameter D[3,2] ( $\mu\text{m}$ ):** 8.993

**CE Diameter STDV ( $\mu\text{m}$ ):** 2.11

**Particles Counted:** 37661

**CE Diameter RSD (%):** 93.66

Volume transformation: CE Diameter ( $\mu\text{m}$ ) smoothed over 11 points



Record 31: Bactiblock 101

## Comments:

Morphologi

Malvern Instruments Worldwide

[www.malvern.com](http://www.malvern.com)