

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Esacol MX155

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Russo

Edited On: 19 March 2014 01:21:27

Date: 19 March 2014 01:21:27

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 10x, 20x

Particles Counted: 1299547

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[n, 0.1] (μm): 0.99

CE Diameter Maximum (μm): 183.62

CE Diameter D[n, 0.5] (μm): 2.44


CE Diameter Mean (μm): 3.87

CE Diameter D[n, 0.9] (μm): 6.53

CE Diameter STDV (μm): 6.24

CE Diameter RSD (%): 161.07

CE Diameter (μm) smoothed over 11 points



Record 308: Esacol MX155

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com