

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: Calce Femia

SOP Name: Calci.vsop

Edited: False

User Name: Domenico

Edited On: 19 December 2014 08:56:15

Date: 19 December 2014 08:56:15

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[v, 0.1]: 5.618

CE Diameter Maximum (μm): 69.59

CE Diameter D[v, 0.5]: 12.36

CE Diameter D[4,3] (μm): 16.15

CE Diameter D[v, 0.9]: 31.68


CE Diameter D[3,2] (μm): 10.46

CE Diameter STDV (μm): 3.34

Particles Counted: 644783

CE Diameter RSD (%): 71.54

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 7: Calce Femia

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com