

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: Anidrite fluorsid

SOP Name: GESSI.vsop

Edited: True

User Name: Antonio Roberto

Edited On: 16 June 2014 13:12:20

Date: 16 June 2014 13:05:50

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[v, 0.1]: 4.323

CE Diameter Maximum (μm): 49.32

CE Diameter D[v, 0.5]: 9.517

CE Diameter D[4,3] (μm): 14.91

CE Diameter D[v, 0.9]: 32.32


CE Diameter D[3,2] (μm): 8.367

CE Diameter STDV (μm): 2.51

Particles Counted: 10124

CE Diameter RSD (%): 64.09

Volume transformation: CE Diameter (μm) smoothed over 11 points



Record 94: Anidrite fluorsid

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com