

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: Eloset 543

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: morfologi

Edited On: 18 April 2014 12:01:45

Date: 18 April 2014 12:01:45

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Particles Counted: 3181

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[n, 0.1] (μm): 0.85

CE Diameter Maximum (μm): 105.44

CE Diameter D[n, 0.5] (μm): 1.97


CE Diameter Mean (μm): 4.79

CE Diameter D[n, 0.9] (μm): 11.95

CE Diameter STDV (μm): 7.99

CE Diameter RSD (%): 166.74

CE Diameter (μm) smoothed over 11 points



Record 19: Eloset 543

Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com