

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: FL 3210

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Morphologi

Edited On: 10 April 2014 14:08:06

Date: 10 April 2014 14:08:06

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 20x

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[v, 0.1]: 20.11

CE Diameter Maximum (μm): 156.51

CE Diameter D[v, 0.5]: 69.76

CE Diameter D[4,3] (μm): 76.22

CE Diameter D[v, 0.9]: 154.5


CE Diameter D[3,2] (μm): 43.12

CE Diameter STDV (μm): 10.21

Particles Counted: 2995

CE Diameter RSD (%): 120.81

Volume transformation: CE Diameter (μm) smoothed over 11 points



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com