

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: Esacol MX155

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Russo

Edited On: 19 March 2014 01:21:27

Date: 19 March 2014 01:21:27

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 10x, 20x

Result

CE Diameter Minimum (μm): 0.54

CE Diameter D[v, 0.1]: 28.73

CE Diameter Maximum (μm): 183.62

CE Diameter D[v, 0.5]: 66.16

CE Diameter D[4,3] (μm): 67.83

CE Diameter D[v, 0.9]: 105.3


CE Diameter D[3,2] (μm): 44.08

CE Diameter STDV (μm): 6.24

Particles Counted: 1299547

CE Diameter RSD (%): 161.07

Volume transformation: CE Diameter (μm) smoothed over 11 points



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com