

# CE Diameter Report - Volume Distribution



## Measurement Details

**Sample Name:** Carbonato di Calcio Femia

**SOP Name:** Calci.vsop

**Edited:** False

**User Name:** Domenico

**Edited On:** 19 December 2014 17:56:38

**Date:** 19 December 2014 17:56:38

## Analysis

**SOP Analysis:** 3.0

**SOP Fill Holes:** True

**SOP Segmentation:** None

**SOP Trash Size:** 10

**SOP Optic(s) used:** 50x

## Result

**CE Diameter Minimum ( $\mu\text{m}$ ):** 0.22

**CE Diameter D[v, 0.1]:** 6.249

**CE Diameter Maximum ( $\mu\text{m}$ ):** 76.58

**CE Diameter D[v, 0.5]:** 15.14

**CE Diameter D[4,3] ( $\mu\text{m}$ ):** 18.53

**CE Diameter D[v, 0.9]:** 35.58


**CE Diameter D[3,2] ( $\mu\text{m}$ ):** 12.06

**CE Diameter STDV ( $\mu\text{m}$ ):** 3.76

**Particles Counted:** 245935

**CE Diameter RSD (%):** 86.50

Volume transformation: CE Diameter ( $\mu\text{m}$ ) smoothed over 11 points



## Comments:

Morphologi

Malvern Instruments Worldwide

[www.malvern.com](http://www.malvern.com)