CE Diameter Report - Number Distribution



Measurement Details

Sample Name: OMYACARB 1T UM --

SOP Name: MATERIE PRIME.vsop Edited: True

User Name: Edited On: 27 May 2014 18:41:37

Date: 27 May 2014 18:34:55

Analysis

SOP Analysis: 3.0 SOP Fill Holes: True

SOP Segmentation: None SOP Trash Size: 10

SOP Optic(s) used: 50x Particles Counted: 21771


Result

CE Diameter Minimum (µm): 0.22 CE Diameter D[n, 0.1] (µm): 1.23

CE Diameter Maximum (µm): 38.08 CE Diameter D[n, 0.5] (µm): 2.96

CE Diameter Mean (µm): 3.51 CE Diameter D[n, 0.9] (µm): 6.43

CE Diameter STDV (µm): 2.28 CE Diameter RSD (%): 65.10



Comments:

Morphologi Malvern Instruments Worldwide www.malvern.com