

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: OMYACARB 1T UM +

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 27 May 2014 12:11:37

Date: 27 May 2014 11:57:01

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Particles Counted: 5821

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[n, 0.1] (μm): 0.42

CE Diameter Maximum (μm): 21.73


CE Diameter D[n, 0.5] (μm): 2.01

CE Diameter Mean (μm): 2.90

CE Diameter D[n, 0.9] (μm): 6.64

CE Diameter STDV (μm): 2.78

CE Diameter RSD (%): 95.81



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com