

CE Diameter Report - Volume Distribution



Measurement Details

Sample Name: FillerGIMOD

SOP Name: MATERIE PRIME.vsop

Edited: False

User Name: Morphologi

Edited On: 28 April 2015 22:20:54

Date: 28 April 2015 22:20:54

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 50x

Result

CE Diameter Minimum (μm): 0.22

CE Diameter D[v, 0.1]: 4.825

CE Diameter Maximum (μm): 36.10

CE Diameter D[v, 0.5]: 11.4

CE Diameter D[4,3] (μm): 12.68

CE Diameter D[v, 0.9]: 22.4


CE Diameter D[3,2] (μm): 8.628

CE Diameter STDV (μm): 2.41

Particles Counted: 184070

CE Diameter RSD (%): 103.40

Volume transformation: CE Diameter (μm) smoothed over 11 points



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com