

CE Diameter Report - Number Distribution



Measurement Details

Sample Name: VINNAPAS 4410 L (Wacker)

SOP Name: MATERIE PRIME.vsop

Edited: True

User Name:

Edited On: 23 May 2014 10:57:16

Date: 23 May 2014 10:50:19

Analysis

SOP Analysis: 3.0

SOP Fill Holes: True

SOP Segmentation: None

SOP Trash Size: 10

SOP Optic(s) used: 10x

Particles Counted: 17637

Result

CE Diameter Minimum (μm): 1.08

CE Diameter D[n, 0.1] (μm): 2.23

CE Diameter Maximum (μm): 187.66


CE Diameter D[n, 0.5] (μm): 7.74

CE Diameter Mean (μm): 11.22

CE Diameter D[n, 0.9] (μm): 22.10

CE Diameter STDV (μm): 12.79

CE Diameter RSD (%): 114.00



Comments:

Morphologi

Malvern Instruments Worldwide

www.malvern.com